

Product Datasheet

SARS Nucleocapsid antibody (orb753269)

Description	Rabbit polyclonal antibody to SARS Nucleocapsid
Reactivity	Virus
Conjugation	Unconjugated
Immunogen	The antibody was developed by immunizing rabbits with synthetic peptides corresponding to amino a.a. 399-411 of putative SARS nucleocapsid (Genbank accession no. YP_009724397.2)
Preservatives	1mg/ml in PBS and 0.02% sodium azide.
Storage	Stability: Store at 4°C, stable for 2 weeks. For long-term storage, store at -20°C
Note	For research use only
Application notes	Applications: Western Blot 1:100-1:2000, Immunohistochemistry, Immunocytochemistry/Immunofluorescence 1:10-1:500
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.