

Product Datasheet

SARS Spike antibody (orb753270)

Description	Rabbit polyclonal antibody to SARS Spike
Reactivity	Virus
Conjugation	Unconjugated
Immunogen	The antibody was developed by immunizing rabbits with synthetic peptides corresponding to amino acids 1124-1140 of the putative SARS Spike glycoprotein (Genbank accession no. NP_828851.1)
Preservatives	0.5mg/ml in PBS, 0.05% rAlbumin & 0.05% sodium azide.
Storage	Stability: Store at 4°C, stable for 2 weeks. For long-term storage, store at -20°C
Note	For research use only
Application notes	Applications: Western Blot 0.05 ug/ml, Immunocytochemistry/Immunofluorescence
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.