



Product Datasheet

ZSCAN4 Antibody (orb75364)

Description	ZSCAN4 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	Q8NAM6
Expiration Date	12 months from date of receipt.



Western blot analysis of mouse lung tiss...



Immunohistochemical analysis of formalin...