



Product Datasheet

ZC3H12B Antibody (orb75392)

Description	ZC3H12B Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	Q5HYM0
Expiration Date	12 months from date of receipt.



Western blot analysis of mouse brain tis...



Immunohistochemical analysis of formalin...