


## Product Datasheet

# Swine H1N1 Nucleocapsid Protein Antibody (orb75405)

<b>Description</b>	Swine H1N1 Nucleocapsid Protein Antibody
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Virus
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA
<b>Storage</b>	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
<b>Note</b>	For research use only
<b>Clonality</b>	Recombinant
<b>Uniprot ID</b>	<a href="#">C4AL25</a>
<b>Expiration Date</b>	12 months from date of receipt.



Bar chart illustrates about the ELISA re...