



Product Datasheet

APC7 Antibody (orb75412)

Description	APC7 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	Q9UJX3
Expiration Date	12 months from date of receipt.



Western blot analysis of rat kidney tiss...



Confocal immunofluorescence analysis of ...