

## Product Datasheet

### Mouse Lymphotactin protein (orb754619)

|                            |                                                                                                                                                  |
|----------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Recombinant mouse Lymphotactin protein                                                                                                           |
| <b>Reactivity</b>          | Mouse                                                                                                                                            |
| <b>Tested Applications</b> | ELISA, WB                                                                                                                                        |
| <b>Preservatives</b>       | Lyophilized from a 0.2 µm filtered solution in 10 mM Hepes, 500 mM NaCl with 5% trehalose, pH 7.4                                                |
| <b>Storage</b>             | -20°C for 12 months as lyophilized; 2-8°C for 1 month under sterile conditions after reconstitution                                              |
| <b>Tag</b>                 | N-terminal His Tag or N-terminal His-IF2DI Tag, determined during production process                                                             |
| <b>Note</b>                | For research use only                                                                                                                            |
| <b>Application notes</b>   | Centrifuge the vial at 10,000 rpm for 1 minute, reconstitute at 200 µg/ml in sterile distilled water by gentle pipetting 2-3 times, don't vortex |
| <b>Purity</b>              | Greater than 95% as determined by SDS-PAGE                                                                                                       |
| <b>Source</b>              | E.Coli                                                                                                                                           |
| <b>MW</b>                  | 10.1 kDa                                                                                                                                         |
| <b>Uniprot ID</b>          | <b>P47993</b>                                                                                                                                    |
| <b>Expiration Date</b>     | 6 months from date of receipt.                                                                                                                   |

---

#### Biorbyt Ltd.

7 Signet Court, Swann's Road,  
Cambridge, CB5 8LA, United Kingdom  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

#### Biorbyt LLC.

68 TW Alexander Drive,  
Durham, NC, 27713, United States  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)