



Product Datasheet

SCRN3 Antibody (orb75462)

Description	SCRN3 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	Q0VDG4
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...



Western blot analysis of 293 cell tissue...