



Product Datasheet

TRIM71 Antibody (orb75475)

Description	TRIM71 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	Q2Q1W2
Expiration Date	12 months from date of receipt.



Western blot analysis of human brain tis...



Immunohistochemical analysis of formalin...