



Product Datasheet

PLAC3 Antibody (orb75511)

Description	PLAC3 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	Q9BXP8
Expiration Date	12 months from date of receipt.



Western blot analysis of HeLa cell lysat...



Immunohistochemical analysis of formalin...