

## Product Datasheet

### RPA32 (phospho-Thr21) antibody (FITC) (orb7552)

<b>Description</b>	Rabbit polyclonal to RPA32 (phospho-Thr21) antibody conjugated to FITC
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Mouse, Rat
<b>Conjugation</b>	FITC
<b>Tested Applications</b>	IF, IHC-Fr, IHC-P, WB
<b>Immunogen</b>	KLH conjugated Synthesised phosphopeptide derived from human RPA2 around the phosphorylation site of Thr21 GY(p-T)Q
<b>Target</b>	RPA2
<b>Preservatives</b>	Lyophilized or Liquid
<b>Form/Appearance</b>	Lyophilized or Liquid
<b>Concentration</b>	- 100 µg (in 100 µl): 1 mg/ml
<b>Storage</b>	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>MW</b>	29 kDa
<b>Uniprot ID</b>	<a href="#">P15927</a>
<b>Dilution Range</b>	IF: 1:50-200
<b>Expiration Date</b>	12 months from date of receipt.