

Product Datasheet

p24 antibody (orb756485)

Description	Mouse monoclonal antibody to p24
Species/Host	Mouse
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant Capsid protein p24
Target	p24
Form/Appearance	liquid
Storage	0.01MPBS pH 7.4,-20°C for 12 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	≥95% as determined by SDS-PAGE
Uniprot ID	P04591
Expiration Date	12 months from date of receipt.