

Product Datasheet

TRIAD3A antibody (orb756626)

Description	Rabbit polyclonal antibody to TRIAD3A
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	E3 ubiquitin-protein ligase RNF216
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3 , -20°C for 12 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	120-130 kDa
Uniprot ID	Q9NWF9
Dilution Range	WB: 1:200-1:2000, IHC-P: 1:20-1:200
Expiration Date	12 months from date of receipt.