

Product Datasheet

2019 nCOV N protein antibody (orb756655)

Description	Mouse monoclonal antibody to 2019 nCOV N protein
Species/Host	Mouse
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, GICA, WB
Immunogen	2019 nCOV N protein
Target	SARS-CoV-2 N protein
Form/Appearance	liquid
Storage	CB buffer, pH 7.5,-20°C for 12 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	≥95% as determined by SDS-PAGE
MW	48 kDa
Dilution Range	ELISA: 0.1-1 μg/ml
Expiration Date	12 months from date of receipt.