

## Product Datasheet

### 2019 nCOV S1 protein antibody (orb756658)

<b>Description</b>	Human polyclonal antibody to 2019 nCOV S1 protein
<b>Species/Host</b>	Human
<b>Reactivity</b>	Virus
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA, WB
<b>Immunogen</b>	S1-RBD
<b>Target</b>	SARS-CoV-2 S1 protein
<b>Form/Appearance</b>	liquid
<b>Storage</b>	PBS pH 7.4, -20°C for 24 months (Avoid repeated freeze / thaw cycles.)
<b>Note</b>	For research use only
<b>Isotype</b>	IgG
<b>Clonality</b>	Monoclonal
<b>Purity</b>	≥95% as determined by SDS-PAGE
<b>Uniprot ID</b>	<a href="#">P0DTC2</a>
<b>Dilution Range</b>	ELISA: 0.1-1 µg/ml, WB: 1:1000-1:10000
<b>Expiration Date</b>	12 months from date of receipt.