


## Product Datasheet

### Human NDUFS1 ELISA Kit (orb756681)

<b>Description</b>	Human NDUFS1 (NADH-ubiquinone oxidoreductase core subunit S1) ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	0.156-10ng/ml
<b>Target</b>	NDUFS1
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Uniprot ID</b>	<b>P28331</b>
<b>Sensitivity</b>	0.094ng/ml
<b>Expiration Date</b>	Please enquire.



Human  
NDUFS1(NADH-  
ubiquinone  
oxidoreduct...