


Product Datasheet

Human SARM1 ELISA Kit (orb756697)

| | |
|------------------------|--|
| Description | Human SARM1 (Sterile alpha and TIR motif-containing protein 1) ELISA Kit |
| Reactivity | Human |
| Range | 0.156-10ng/ml |
| Target | SARM1 |
| Storage | 2-8 °C for 6 months |
| Note | For research use only |
| Uniprot ID | Q6SZW1 |
| Sensitivity | 0.094ng/ml |
| Expiration Date | Please enquire. |



Standard Curve of Human SARM1 ELISA Kit.