


## Product Datasheet

### Human TSG101 ELISA Kit (orb756714)

<b>Description</b>	Human TSG101 (Tumor susceptibility gene 101 protein) ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	31.25-2000pg/ml
<b>Target</b>	TSG101
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Uniprot ID</b>	<b>Q99816</b>
<b>Sensitivity</b>	18.75pg/ml
<b>Expiration Date</b>	Please enquire.



Human  
TSG101  
(Tumor  
susceptibility  
gene ...