


Product Datasheet

Human HAVCR2 ELISA Kit (orb756728)

Description	Human HAVCR2 (Hepatitis A virus cellular receptor 2) ELISA Kit
Reactivity	Human
Range	0.156-10ng/ml
Target	HAVCR2
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	Q8TDQ0
Sensitivity	0.094ng/ml
Expiration Date	Please enquire.



Human
HAVCR2(Hepatitis
A virus cellular ...