


## Product Datasheet

### Human CDCP1 ELISA Kit (orb756767)

<b>Description</b>	Human CDCP1 (CUB domain-containing protein 1) ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	0.156-10ng/ml
<b>Target</b>	CDCP1
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Uniprot ID</b>	<b>Q9H5V8</b>
<b>Sensitivity</b>	0.094ng/ml
<b>Expiration Date</b>	Please enquire.



Human  
CDCP1(CUB  
domain-  
containing  
protei...