


## Product Datasheet

### Human CLEC4M ELISA Kit (orb756769)

<b>Description</b>	Human CLEC4M (C-type lectin domain family 4 member M) ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	3.125-200ng/ml
<b>Target</b>	CLEC4M
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Uniprot ID</b>	<b>Q9H2X3</b>
<b>Sensitivity</b>	1.875ng/ml
<b>Expiration Date</b>	Please enquire.



Human CLEC4M(C-type lectin domain family...