


Product Datasheet

Rabbit Alb ELISA Kit (orb756797)

Description	Mouse Alb (albumin) elisa kit
Reactivity	Mouse
Range	7.813-500ng/ml
Target	Alb
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Uniprot ID	P07724
Sensitivity	4.688ng/ml
Expiration Date	Please enquire.



Mouse
Alb(albumin)
elisa kit