


Product Datasheet

LPS ELISA Kit (orb756819)

Description	LPS (Lipopolysaccharides) ELISA Kit
Reactivity	All
Range	0.313-20ug/ml
Target	LPS
Storage	2 - 8 °C for 12 months
Note	For research use only
Application notes	Competitive ELISA, Coated with Antibody
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Sensitivity	0.188ug/ml
Expiration Date	Please enquire.



LPS(Lipopolysaccharides) ELISA Kit

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)