



Product Datasheet

Ly6G antibody (Biotin) (orb756847)

Description	Rat monoclonal antibody against Ly6G conjugated to Biotin
Reactivity	Mouse
Conjugation	Biotin
Tested Applications	FC, IHC-Fr, IHC-P, IP, WB
Immunogen	Murine granulocytes
Target	Ly6G
Preservatives	Stabilizing phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: Recommended dilution: 2-5 µg/ml.
Isotype	Rat IgG2b
Clonality	Monoclonal
Purity	Purified antibody is conjugated with biotin LC-NHS ester under optimum conditions and unconjugated antibody and free biotin are removed by size-exclusion chromatography.
Uniprot ID	A0A087WQF5
Entrez	546644
Dilution Range	Flow cytometry: Recommended dilution: 2-5 µg/ml.
Expiration Date	12 months from date of receipt.



Separation of murine myeloid Ly6G posi...



Flow cytometry surface staining pattern ...