
Product Datasheet

Human Galectin 9 Protein, hFc Tag (orb757445)

Description

Recombinant human Galectin 9 Protein

| | |
|------------------------|---|
| Reactivity | Human |
| Conjugation | Unconjugated |
| Target | Galectin 9 |
| Preservatives | Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions. |
| Form/Appearance | Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution. |
| Storage | Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature. |
| Tag | N-Human Fc Tag |
| Note | For research use only |
| Purity | The purity of the protein is greater than 95% as determined by SDS-PAGE and Coomassie blue staining. |
| Source | Mammalian |
| MW | The protein has a predicted molecular mass of 61.9 kDa after removal of the signal peptide. The apparent molecular mass of hFc-Galectin 9 is approximately 70-100 kDa due to glycosylation. |
| Uniprot ID | 000182 |
| Expiration Date | 6 months from date of receipt. |