

## Product Datasheet

### PGN9856 (orb757618)

<b>Description</b>	PGN-9856 (PGN9856) is a novel potent, selective prostaglandin EP2 receptor agonist with pKi of > 8.3...
<b>Note</b>	For research use only
<b>Formula</b>	$C_{21}H_{15}F_2NO_4$
<b>MW</b>	383.351
<b>CAS Number</b>	[863704-91-0]
<b>SMILES</b>	<chem>C(O)(=O)COC1=CC=CC(NC(C2C(F)=CC=C(C3=CC=CC(F)=C3)C=2)=O)=C1</chem>
<b>Expiration Date</b>	12 months from date of receipt.