

Product Datasheet

TAO Kinase inhibitor 1 (orb757621)

Description	TAOK inhibitor 43 is a novel potent, selective, ATP-competitive TAO Kinase (TAOK) inhibitor with IC5...
Note	For research use only
Formula	$C_{25}H_{24}N_2O_2$
MW	384.47
CAS Number	[850467-66-2]
SMILES	<chem>O=C(CNC(C1=CC=C(C2=CC=CC=C2)C=C1)=O)NC(CCC3)C4=C3C=CC=C4</chem>
Expiration Date	12 months from date of receipt.