

Product Datasheet

3'3'-cGAMP (sodium salt) (orb757686)

Description	3'3'-cGAMP Fluorinated (c-[2' ³ FdGMP]-[2' ³ FdAMP]) is a synthetic analog of cyclic guanosine monopho...
Note	For research use only
Formula	C ₂₀ H ₂₄ N ₁₀ O ₁₃ P ₂
MW	674.4112844
CAS Number	[849214-04-6]
SMILES	<chem>O[C@@H]1[C@@H]2OP(OC[C@H]3O[C@@H](N4C=NC5C(=NC=NC4=5)N)[C@H](O)[C@@H]3OP(OC[C@H]2O[C@H]1N1C=NC2C(N=C(NC1=2)N)=O)(O)=O)(O)=O</chem>
Expiration Date	12 months from date of receipt.