

Product Datasheet

CD73 inhibitor AB-680 (orb757691)

Description	AB680 (AB-680) is a highly potent, selective, reversible inhibitor of CD73 with Ki/IC50 of 4.9/70 pM...
Note	For research use only
Formula	$C_{20}H_{24}ClFN_4O_9P_2$
MW	580.8248892
CAS Number	[2105904-82-1]
SMILES	<chem>C1C1C([H])=C(C2C([H])=NN(C=2N=1)[C@@]1([H])[C@@]([H])([C@@]([H])([C@@]([H])(C([H])([H])OP(C([H])([H])P(=O)(O[H])O[H])(=O)O[H])O1)O[H])O[H])N([H])[C@@]([H])(C([H])([H])[H])C1=C([H])C([H])=C([H])C([H])=C1F</chem>
Expiration Date	12 months from date of receipt.