



Product Datasheet

ABCA7 Antibody (orb75805)


Description	ABCA7 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	Q8IZY2
Expiration Date	12 months from date of receipt.



Western blot analysis of 293 cell lysate...



Immunohistochemical of human spleen tiss...



Immunofluorescence analysis of 293 cells...