

Product Datasheet

EFNB1/2 (Phospho-Tyr330) antibody (orb758122)

Description	Rabbit polyclonal antibody to EFNB1/2 (Phospho-Tyr330)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	C-terminal
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P98172, P52799
Expiration Date	12 months from date of receipt.