

Product Datasheet

c-Kit (Phospho-Tyr703) antibody (orb758140)

Description	Rabbit polyclonal antibody to c-Kit (Phospho-Tyr703)
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Internal
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P10721
Expiration Date	12 months from date of receipt.