

Product Datasheet

14-3-3 Zeta/Delta (Phospho-Thr232) antibody (orb758154)

Description	Rabbit polyclonal antibody to 14-3-3 Zeta/Delta (Phospho-Thr232)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	C-terminal
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P63104
Expiration Date	12 months from date of receipt.