

Product Datasheet

PKC-pan (Phospho-Thr497) antibody (orb758159)

| | |
|----------------------------|---|
| Description | Rabbit polyclonal antibody to PKC-pan (Phospho-Thr497) |
| Species/Host | Rabbit |
| Reactivity | Human, Mouse, Rat |
| Conjugation | Unconjugated |
| Tested Applications | IF, IHC, WB |
| Immunogen | C-terminal |
| Concentration | 1mg/ml |
| Storage | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| Note | For research use only |
| Clonality | Polyclonal |
| Uniprot ID | P17252, P05129, P05771, Q02156, P24723, Q05655, Q04759, Q05513 |
| Expiration Date | 12 months from date of receipt. |