

Product Datasheet

Catenin- beta (Phospho-Tyr654) antibody (orb758168)

Description	Rabbit polyclonal antibody to Catenin-beta (Phospho-Tyr654)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	C-terminal
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P35222
Expiration Date	12 months from date of receipt.