

Product Datasheet

p57 Kip2 (Phospho-Thr310) antibody (orb758183)

Description	Rabbit polyclonal antibody to p57 Kip2 (Phospho-Thr310)
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	C-terminal
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P49918
Expiration Date	12 months from date of receipt.