

Product Datasheet

CDC16/APC6 (Phospho-Ser560) antibody (orb758188)

Description	Rabbit polyclonal antibody to CDC16/APC6 (Phospho-Ser560)
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	C-terminal
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Uniprot ID	Q13042
Expiration Date	12 months from date of receipt.