

## Product Datasheet

# GIRK1/KIR3.1/KCNJ3 (Phospho-Ser185) antibody (orb758189)

<b>Description</b>	Rabbit polyclonal antibody to GIRK1/KIR3.1/KCNJ3 (Phospho-Ser185)
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human, Mouse, Rat
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	IF, IHC, WB
<b>Immunogen</b>	Internal
<b>Concentration</b>	1mg/ml
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Clonality</b>	Polyclonal
<b>Uniprot ID</b>	<a href="#">P48549</a>
<b>Expiration Date</b>	12 months from date of receipt.