

## Product Datasheet

# IL-8R beta/CDw128 beta (Phospho-Ser347) antibody (orb758198)

<b>Description</b>	Rabbit polyclonal antibody to IL-8R beta/CDw128 beta (Phospho-Ser347)
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human, Mouse
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	IF, IHC, WB
<b>Immunogen</b>	C-terminal
<b>Concentration</b>	1mg/ml
<b>Note</b>	For research use only
<b>Clonality</b>	Polyclonal
<b>Uniprot ID</b>	<b>P25025</b>
<b>Expiration Date</b>	12 months from date of receipt.