

Product Datasheet

AurB/C (Phospho-Thr236/202) antibody (orb758216)

Description	Rabbit polyclonal antibody to AurB/C (Phospho-Thr236/202)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	Internal
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Uniprot ID	Q9UQB9 , Q96GD4
Expiration Date	12 months from date of receipt.