

Product Datasheet

Telomerase (Phospho-Ser824) antibody (orb758224)

Description	Rabbit polyclonal antibody to Telomerase (Phospho-Ser824)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Internal
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Uniprot ID	O14746
Expiration Date	12 months from date of receipt.