

Product Datasheet

APC1 (Phospho-Ser355) antibody (orb758262)

Description	Rabbit polyclonal antibody to APC1 (Phospho-Ser355)
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	Internal
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Uniprot ID	Q9H1A4
Expiration Date	12 months from date of receipt.