


Product Datasheet

Hemagglutinin Monoclonal Antibody (orb75864)

Description	Hemagglutinin Monoclonal Antibody
Species/Host	Mouse
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Monoclonal
Uniprot ID	Q692M2
Expiration Date	12 months from date of receipt.



Western blot analysis of 25 ng of recomb...