

Product Datasheet

IgM antibody (orb758769)

Description	Rabbit monoclonal antibody to IgM
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	Blocking, ELISA, FC, IHC-Fr, IP
Immunogen	BALB/c mice immunized with human IgM.
Target	IgM
Preservatives	PBS with 0.02% Proclin 300.
Concentration	1 mg/ml
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Application notes	This antibody binds to IgM, which acts in primary immune defines and is involved in early recognition of antigens.
Isotype	IgG
Clonality	Monoclonal
Purity	Purified
Uniprot ID	P01871
Expiration Date	12 months from date of receipt.