


Product Datasheet

Seasonal H1N1 Neuraminidase Monoclonal Antibody (orb75890)

Description	Seasonal H1N1 Neuraminidase Monoclonal Antibody
Species/Host	Mouse
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Monoclonal
Uniprot ID	B1AGS8
Expiration Date	12 months from date of receipt.



Seasonal influenza A N1 antibody specifi...