

Product Datasheet

Beta 2 adrenergic receptor antibody (orb758965)

Description	Rabbit monoclonal antibody to Beta 2 adrenergic receptor
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IF, IP, WB
Immunogen	This antibody was raised by immunizing mice with beta-2 adrenergic receptor.
Target	Beta 2 adrenergic receptor
Preservatives	PBS with 0.02% Proclin 300.
Concentration	1 mg/ml
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Application notes	This antibody is recommended for the detection and analysis of the beta-2 adrenergic receptor. It was successfully used in the in vivo phosphorylation assays during Western blot detections of beta-2 adrenergic receptor (Huang et al., 2009).
Isotype	IgG
Clonality	Monoclonal
Purity	Purified
Uniprot ID	P07550
Expiration Date	12 months from date of receipt.