
Product Datasheet

Beta 2 adrenergic receptor antibody (orb758967)

Description

Rabbit monoclonal antibody to Beta 2

Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IP, WB
Immunogen	This antibody was raised by immunizing mice with beta-2 adrenergic receptor.
Target	Beta 2 adrenergic receptor
Preservatives	PBS with 0.02% Proclin 300.
Concentration	1 mg/ml
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Application notes	This antibody is recommended for the identification and analysis of the β 2-adrenergic receptor. For instance, it was used to immunostain vascular smooth muscle cells (VSMCs) against beta-2-adrenergic receptor to determine the fate of beta-2ARs ubiquitinated in a carvedilol responsive manner (Han et al., 2012). Furthermore, the antibody was utilized in the quantification of the surface level of β -2ARs via ELISA (Han et al., 2012).
Isotype	IgG
Clonality	Monoclonal
Purity	Purified
Uniprot ID	P07550
Expiration Date	12 months from date of receipt.