

Product Datasheet

Primate TSLP protein (orb759100)

Description	Recombinant Cynomolgus Thymic Stromal Lymphopoietin is produced by our E.coli expression system and ...
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Form/Appearance	Lyophilized from a 0.2 um filtered solution of 20mM PB , 500mM NaCl, pH 7.4.
Storage	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
Tag	C-6His
Note	For research use only
Uniprot ID	XP_005557555.1
Expiration Date	6 months from date of receipt.