

Product Datasheet

Mouse Siglec-15 protein (orb759153)

Description	Recombinant Mouse Sialic Acid-binding Ig-like lectin 15 is produced by our Mammalian expression syst...
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Form/Appearance	Lyophilized from a 0.2 um filtered solution of PBS,150mMnacl,5%Thehalose,0.3%chaps,Ph7.4.
Storage	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
Tag	C-6His
Note	For research use only
Uniprot ID	A7E1W8
Expiration Date	6 months from date of receipt.