

Product Datasheet

Human IL23A & IL12B Heterodimer protein (orb759162)

Description	Recombinant Chimeric Human&Mouse Interleukin-23 is produced by our Mammalian expression system and t...
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Form/Appearance	Lyophilized from a 0.2 um filtered solution of PBS,pH7.4.
Storage	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
Tag	C-6His
Note	For research use only
Uniprot ID	Q9NPF7, P43432
Expiration Date	6 months from date of receipt.